

# Reichert® Phoroceptor® VRx & SightChek™ Wireless Dongle **16233**

---



**Reichert, Inc.**

www.reichert.com • 3362 Walden Avenue, Suite 100, Depew, NY 14043 USA

Toll Free: 888-849-8955 • Tel: +1 716-686-4500 • Fax: +1 716-686-4545

**16233-101 Rev. C** (November 26, 2018)

©2018 AMETEK, Inc & Reichert, Inc. Bluetooth trademark is owned by Bluetooth SIG. **CE**



Welcome to the Reichert® Phoroptor® VRx **Wireless Dongle 16233**, allowing Bluetooth® wireless connectivity between Phoroptor® VRx and numerous ophthalmic devices. To use the Wireless Dongle with another device, an additional Reichert adapter kit is required. Please contact Reichert or your authorized Reichert distributor to determine the compatible kit for your device.

**NOTE:** This document will familiarize you with your new Phoroptor VRx accessory. *To use this Wireless Dongle you must first configure it to communicate with the Phoroptor VRx. Please refer to the Phoroptor VRx Digital Refraction System User's Guide (16241-101) for complete installation, configuration, and compliance information.* To obtain a copy of 16241-101, please contact Reichert, your Reichert authorized dealer, or you may download it from the Phoroptor VRx product page at [www.reichert.com](http://www.reichert.com). Actual components may differ in appearance from those depicted here.



## Wireless Dongle Contains:

---

- 1 Wireless Dongle AC-Power Supply w/ Plug Adapter
- 2 Wireless Dongle
- 3 Wireless Dongle DB9 Programming Cable