

# Reichert® Phoroceptor® VRx Wireless Dongle 16233

---



Reichert, Inc.  
www.reichert.com • 3362 Walden Avenue, Depew, NY 14043 USA  
Toll Free: 888-849-8955 • Tel: +1 716-686-4500 • Fax: +1 716-686-4545  
**16233-101 Rev. A** (December 17, 2014)  
©2014 AMETEK, Inc & Reichert, Inc. Bluetooth trademark is owned by Bluetooth SIG. CE



Welcome to the Reichert® Phoroceptor® VRx **Wireless Dongle 16233**, allowing Bluetooth® wireless connectivity between **Phoroceptor® VRx** and numerous ophthalmic devices. To use the Wireless Dongle with another device, an additional Reichert adapter kit is required. Please contact Reichert or your authorized Reichert distributor to determine the compatible kit for your device.

**NOTE:** This document will familiarize you with your new Phoroceptor VRx accessory. *To use this Wireless Dongle you must first configure it to communicate with the Phoroceptor VRx. Please refer to the Phoroceptor VRx Digital Refraction System User's Guide (16241-101) for complete installation, configuration, and compliance information.* To obtain a copy of 16241-101, please contact Reichert, your Reichert authorized dealer, or you may download it from the Phoroceptor VRx product page at [www.reichert.com](http://www.reichert.com). Actual components may differ in appearance from those depicted here.



## Wireless Dongle Contains:

---

- 1 Wireless Dongle
- 2 Wireless Dongle AC-Power Supply w/ Plug Adapter
- 3 Wireless Dongle DB9 Programming Cable